

ISO9001,ISO13485 Quality Management

Diatek

DRW-320 CE **Microplate Washer**



DESIGN BY
DIATEK



Wuxi Hiwell Diatek Instruments Co.,Ltd

Main Features

- ◆ Large LCD display, menu designed for convenient use
- ◆ Whole plate washing or single strip washing
- ◆ 12-way and 8-way manifold included
- ◆ Low residual volume by two pipettes
- ◆ Complete bottom washing
- ◆ Shaking and dipping function
- ◆ Automatic monitoring of vacuum and pressure ,automatic rinse cycle
- ◆ Flat, V-bottom or U-bottom plate and strips washing
- ◆ Large memory to store up to 120 user programmed wash protocols
- ◆ Emergency stop available and liquid warning



Technical Specification

Input	8 key membrane keyboard
Display	High brightness LCD display screen
Manifold	8 pins and 12 pins
Wash mode	Strip mode and plate mode
Wash strips	1~12 adjustable
Applicable well	Flat, V-bottom or U-bottom
Washing times	0~99 times adjustable
Dipping time	0~3600s adjustable
Shaking time	0~600s adjustable
Residual volume	≤1μL/well
Liquid volume	50~3000ul/well
Sipping time	0.1~9.9s
Storage capacity	More than 120 user defined wash board procedures
Power supply	100~240V, 50~60Hz
Dimensions	448(L)×382(W)×163(H)mm
Net weight	8kg
Operating environment	Temperature 5℃~40℃, Max humidity 80%



Wuxi Hiwell Diatek Instruments Co.,Ltd

Tel:+86 510 85160879
Fax:+86 510 85166536

E-mail:sales@diateklab.com
Http:www.diateklab.com

Add:Room 1901~1904,19F,Building2,2166
Taihu West Avenue,Wuxi,China

